

BURN FREE BURNS DRESSING 20 X 20CM



Product Code

CM0333

Available in Sizes

20X20CM

Available in Colours



White
CM0333

Click Medical have formed a partnership with BurnFree, one of the world's leading manufacturers of an innovators in Hydrogel technology and burns first aid products.

BurnFree works!

How? BurnFree gel draws the heat out of burn injuries and into the hydrogel. This cools the burn injury providing pain relief and helps stop the burn progression. BurnFree also keeps the wound site moist so that the body will not lose more water, helping to keep burn victims from entering into shock. BurnFree is a viscous gel. This allows the exposed nerve endings to be covered by the gel, offering near immediate pain relief.

While water will cool the burn, it is often difficult if not impossible to keep water over the entire burn area. Also, a large volume of water is required to cool even a small burn. BurnFree works best because the BurnFree gel is viscous, meaning that it will stay on the burn and draw out the heat, making it easier to cool the burn with BurnFree than with water. By covering the burn site, BurnFree will also help to prevent contamination until medical care is obtained.

Dressings - All Burnfree dressings are designed for first aid use on all types of burn injuries. The dressings are constructed of a special soft, open cell foam and are elastic enough to allow possible swelling or the need to flex joints. Burnfree dressings are non-adherent and the inert foam will not shed threads, strands or fibres that could get embedded in the wound. The pliable foam holds 50 times its weight in gel, meaning the gel will always stay in contact with the burn.

