

## PLASTIC PERFORATED TAPE 2.5CM X 10M

Product Code

CM0430

Available in Sizes

**ONE SIZE**

- Size: 2.5 cm x 9.1 m
- Non Sterile
- Premier Claripore is a transparent Polyethylene based tape material with acrylic based adhesive.
- Perforated and tears easily in clean straight lines, cross-wise and length-wise without fraying edges.
- Ideal for use where tape width must be customised.
- Breathable, making it gentle on the skin and reducing the risk of maceration.
- The tape is easy to remove and leaves minimal adhesive residue.
- Can be written on, sticks well to itself and easy to use with gloves on.
- Always ensure skin is clean and dry before applying tapes, clip hair on the tape adhesion skin area if necessary. Always select a breathable tape where repeated taping is required, or the skin is fragile, this will help to reduce trauma/damage to the skin.
- Consider using a stretchable tape if there is a risk of swelling.
- Do not encircle an entire limb with tape, subsequent swelling may affect circulation. Replace or reposition tape if swelling, oedema or loss of tension in dressing occurs.
- Intended Uses; Light strapping e.g immobilising fingers and toes. Securing dressings. Anchoring tubes and devices (catheters, probes cannulae and nasogastric).
- Latex Free.
- Single Use.
- Store in a cool dry place out of direct sunlight.

