



**BEESWIFT**  
FOCUSED ON SAFETY

Product Code  
CM1650

Available in Sizes  
**ONE SIZE**

Available in Colours  
 White  
CM1650

# DATASHEET


## ZOLL AED 3 CPR UNI-PADZ




- For use exclusively with the ZOLL AED 3®, CPR Uni-padz electrodes have an industry-leading shelf life of five years and can be used on both adult and child patients, so there is no need to maintain two separate sets of pads.
- The ZOLL AED 3 full-colour LCD display includes graphics that coincide with the electrode shapes and colours to support proper pad placement on both adults and children.
- By simply repositioning the pads and activating the Child Mode on the ZOLL AED 3, CPR Uni-padz can be used to treat child patients 8 years of age or younger or under 25 kg (55 lbs.).

### CPR Uni-padz™


1



2



3



OPERATING TEMPERATURE: 0°C to 50°C (32°F to 122°F)  
SHORT TERM STORAGE TEMPERATURE: -30°C to 65°C (-22°F to 149°F)  
LONG TERM STORAGE TEMPERATURE: 0°C to 35°C (32°F to 95°F)

Rx ONLY 32°F 39°C  
NOT MADE WITH NATURAL RUBBER LATEX  
CE 0123

R2138-30 Rev. B Patent: www.zoll.com/patents

**ZOLL**  
ZOLL Medical Corporation  
2400 AED Road  
Chatham, MA 01924-1001 USA  
800.348.9011  
Made in USA

ZOLL International Holding B.V.  
Neukirchweg 18 6662 PV EIJH  
The Netherlands Tel: +31 (0) 481 366410

LOT 2321  
REF 8900-000270  
2026-09-05

(01)00847946024161(17)260905(10)2321

