

INSTANT WARM PERINEAL COMPRESS 5"X 16"



Product Code

RA94640

Available in Sizes

16"

Available in Colours

White
RA94640

- Direct To Skin
- Fold-to-Pop Mechanism (FTP)
- Instant Activation
- One Time Use
- Warm

The Rapid Relief Instant Warm Perineal Compress (Direct-to-skin) quickly delivers soothing warmth, helping to encourage relaxation of tissues and increased blood circulation in the perineal area prior to childbirth, providing maximum patient comfort with minimal nursing intervention.

Features

- Temperature: Delivers approximately 45.6'C/114'F of warm therapy when activated at room temperature
- Extra absorbent compress comes with an adhesive strip to keep it firmly in place
- Longer front-to-back area for better coverage and protection from leakage
- Much more convenient than a warm sitz bath
- FTP (fold-to-pop) technology, allowing activation of the warm compress simply by folding it in half
- Direct-to-skin application

Application

- **INSTANT WARM THERAPY RELIEF:**
- Encourage relaxation of tissues prior to childbirth
- Start the healing process after childbirth and vasectomy procedures (usually 72 hours after the use of cold therapy)
- Muscle tension and menstrual cramps

